



HORSEMAN
of Royal Berkshire



CLASS 1 MOBILITY SCOOTERS

Class 1 mobility scooters are compact and lightweight, easy to transport, usually dismantle for ease of transportation, and are easy to stow. This makes them perfect for taking on holiday, days away or shorter trips at home.

The average weight capacity of Class 1 mobility scooter is around 127 kilograms and will travel around 10 miles on a single charge at a maximum speed of 4mph.



HORSEMAN
of Royal Berkshire



CLASS 2 MOBILITY SCOOTERS

Class 2 mobility scooters are generally more stable and sturdier than the [class one scooter](#) and are able to carry more weight. Although they can also travel further distances, they are, however, generally not as simple to dismantle, are heavier and more awkward to transport.

The maximum user weights are around 133 kilograms to 158 kilograms and can travel speeds of 4mph to 6mph.



HORSEMAN
of Royal Berkshire



ELECTRIC WHEELCHAIRS

An electric powered wheelchair, rather than a mobility scooter, may be a more suitable option, especially if there is a need to stay in the wheelchair throughout the day, or have restricted mobility and require assistance around the home, or are unable to transfer on and off a seat easily.

Some of the smaller electric wheelchairs will dismantle to fit into the smallest of car boots. However, the larger, heavier chairs do not fold or dismantle and are too heavy to lift.



HORSEMAN
of Royal Berkshire



CARER PROPELLED WHEELCHAIRS

Transit wheelchairs are designed to be pushed by another person.

With a smaller rear wheel, Transit wheelchairs are easier to manoeuvre and more compact. Most of our transit or attendant propelled wheelchairs have a simple folding mechanism so that the sides close in.



HORSEMAN
of Royal Berkshire



SELF PROPELLED WHEELCHAIRS

Self-Propelled wheelchairs are designed with a larger wheel at the back to allow the user to push themselves.

Self-propelled wheelchairs give flexibility – they can be pushed by the wheelchair user, but they also easily fold down if you need to stow them away.